

What is most important to me
when closing wounds?



Superior Flexibility and Optimal Strength.

www.unimed.com.sa

REF SS-035T
SurgiSeal™
TOPICAL SKIN ADHESIVE
2-Octyl Cyanoacrylate

SurgiSeal™

TOPICAL SKIN ADHESIVE
2-Octyl Cyanoacrylate

Superior Flexibility and Optimal Strength

SurgiSeal is a topical skin adhesive that provides the optimal balance between strength and flexibility and can be used to close lacerations, wounds, and surgical incisions. SurgiSeal offers all of the benefits of a traditional topical tissue adhesive, as well as:

High-strength 2 Octyl-Cyanoacrylate formulation

Octyl-cyanoacrylate has the greatest clinical tensile strength among available tissue adhesives and a broader indication for clinical use.

High Moisture Vapor Transmission Rate (MVTR)

The proprietary formulation is designed to provide optimal permeability, increasing moisture and oxygen transmission to the wound.

OctylFlex™ Technology

Provides increased flexibility designed to contribute to patient comfort and help preserve the integrity of the adhesive, preventing premature sloughing.

Precision Applicator

SurgiSeal's patented applicator is designed to provide greater coverage with less adhesive material and requires no catalyst to activate its adhesive properties. This self-contained applicator is designed to increase efficiency in the OR, to be used on multiple lacerations, and to reduce waste.

| Description | Product Code | Fill Volume (ml) | Packaging |
|----------------------------------|--------------|------------------|------------------------|
| SurgiSeal™ Topical Skin Adhesive | SS-035T | 0.35ml | 12 Application per box |

UNIMED

United Medical Industries Co. Ltd.

Head Office P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel : +966 1 265 0665
Fax: +966 1 265 0569
www.unimed.com.sa
general@unimed.com.sa

Higher Permeability Promotes Faster Healing

Studies demonstrate that adhesives with high Moisture Vapor Transmission Rate (MVTR) contribute to improved wound care, and tests show that SurgiSeal's MVTR is more than twice that of its leading competitor.

Comparison of Permeability

