

**Yes!**  
It is polyglactin 910

**UNIGLACTIN**

**3 Weeks**  
Effective Tensile  
Strength

**UNIMED<sup>®</sup>**  
SUTURES


United Medical Industries Co. Ltd.

[www.unimed.com.sa](http://www.unimed.com.sa)

# UNIGLACTIN

- Coated Braided Polyglactin 910 B.P Violet and Undyed.
- **Synthetic absorbable sterile surgical suture composed of a copolymer made from 90% glycolide and 10% L-lactide.**
- Available in range of gauge sizes and lengths, non-needled or attached to stainless steel needles of varying types and sizes.
- Minimal tissue reaction.
- Easy passage through tissue.
- Secure Knot tying and easy knot tie down.
- Intended for use where an absorbable suture or ligature is indicated. It can be used in general soft tissue approximation or ligation and in ophthalmic procedures.
- Sterilised by Ethylene Oxide.
- Coating – Coated with mixture of copolymer of glycolide and lactide and Calcium Stearate.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United State Pharmacopoeia (U.S.P.) except for an occasional slight oversize in some gauges.

**UNIGLACTIN - Approximate Breaking Strength Retention**



## Head Office

P.O.Box 51871  
Riyadh 11553 Saudi Arabia  
Tel. +966 1 2650665  
Fax. +966 1 2650569  
general@unimed.com.sa

[www.unimed.com.sa](http://www.unimed.com.sa)