



UNIPRO

UNIMED[®]
SUTURES

United Medical Industries Co. Ltd.
www.unimed.com.sa

UNIPRO

- Monofilament Polypropylene suture B.P. Blue.
- Synthetic non-absorbable sterile surgical suture.
- Available in range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- Inert – Minimal tissue reaction.
- Very smooth for easy passage through tissue.
- Less memory, very pliable- good handling.
- Easy knot tying, excellent knot security.
- Even in infected tissue, polypropylene provides long term strength retention.
- Intended for use in procedures where a smooth monofilament suture offering very long term support is required. It can be used in general soft tissue approximation or ligation. Also, recommended for cardiovascular, ophthalmic and neurological procedures.
- Sterilised by Ethylene Oxide.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), the European Pharmacopoeia (Ph.Eur.), and the United State Pharmacopoeia (U.S.P.).

Head Office

P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel. +966 1 2650665
Fax. +966 1 2650569
general@unimed.com.sa

www.unimed.com.sa