











UNISTEEL MONOFILAMENT STAINLESS STEEL WIRE SUTURE B.P.

NEEDLE		LENGTH	GAUGE
TEMPORARY PACEMAKER WIRES - STANDARD			32 swg
16mm			
Curved round bodied 		45cm	DW 1904D
63mm			
Straight cutting 		if break point required	B1904D
25mm			
1/2 circle round bodied 		45cm	DW 1920
25mm			
1/2 circle round bodied 		45cm	DW 1924D
90mm			
Curved cutting 		if break point required	B1924D

32 swg soft braided, specially insulated stainless steel wire. 7.5cm uninsulated at either end. 30cm of insulation. 45cm overall length. The 63mm needle used can be supplied with a 'breakpoint' if required. 'Breakpoint' needle pacemaker wires are supplied on 59.5cm wire (52 cm of insulation and 7.5cm uninsulated at the non-'breakpoint' end).

UNISTEEL MONOFILAMENT STAINLESS STEEL WIRE SUTURE B.P.

NEEDLE		LENGTH	GAUGE	
MONOFILAMENT STEEL WIRE SUTURE			1 (4)	5 (7)
40mm				
1/2 circle taper cut - heavy Multipack (4 sutures per envelope) 		45cm	DT 1714M	
55mm				
1/2 circle Trocar point - heavy 		75cm	DW 1798H	