



UNICRYL M

6 Weeks
Effective Tensile
Strength

UNIMED[®]
SUTURES


United Medical Industries Co. Ltd.

www.unimed.com.sa

UNICRYL M

- Monofilament Polydioxanone E.P. Violet and Undyed.
- Synthetic absorbable sterile surgical suture prepared from polyester p-dioxanone.
- Available in range of gauge sizes and lengths, non-needled or attached to stainless steel needles of varying types and sizes.
- Minimal tissue reaction.
- Easy passage through tissue.
- Intended for use where an absorbable suture or ligature is indicated. It can be used in all types of tissue approximation, pediatric cardiovascular tissue and ophthalmic surgery.
- Sterilised by Ethylene Oxide.
- Sutures are guaranteed to meet the requirements of the British Pharmacopoeia (B.P.), European Pharmacopoeia (Ph.Eur.), and the United States Pharmacopoeia (U.S.P.).

UNICRYL M - Approximate Breaking Strength Retention



Head Office

P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel. +966 1 2650665
Fax. +966 1 2650569
general@unimed.com.sa