



UNICAPRONE

7 Days
Effective Tensile
Strength

UNIMED[®]
SUTURES


United Medical Industries Co. Ltd.

www.unimed.com.sa

UNICAPRONE

- Monofilament Polyglecaprone 25 B.P. Violet and Undyed.
- Synthetic absorbable sterile surgical suture composed of Poly (glycolide-co-ε-caprolactone).
- Available in range of gauge sizes and lengths, non-needed or attached to stainless steel needles of varying types and sizes.
- Minimal tissue reaction.
- Easy passage through tissue.
- Intended for use where an absorbable suture or ligature is indicated for short term support. It can be used in general soft tissue approximation or ligation.
- Sterilised by Ethylene Oxide.
- Sutures are guaranteed to meet the requirements of British Pharmacopoeia (B.P.) the European Pharmacopoeia (Ph.Eur.), and the United States Pharmacopoeia (U.S.P).

UNICAPRONE - Approximate Breaking Strength Retention



Head Office

P.O.Box 51871
Riyadh 11553 Saudi Arabia
Tel. +966 1 2650665
Fax. +966 1 2650569
general@unimed.com.sa

www.unimed.com.sa